

scanX

ACCESSORIES FOR ALLPRO SCANX SYSTEMS

ALLPRO Imaging offers convenient and efficient accessories specifically designed for your ScanX system. Our durable carrying options make transporting the ScanX a breeze. And, innovative Wall Mounts save precious space in your practice.

SOFT CARRY CASE

Transport and protect your valuable investment with our durable, lightweight Soft Carry Case featuring a contemporary contoured design. Strong, easy-grip handles and a padded, comfortable shoulder strap make it easy to lift and carry. Thick, foam walls ensure extra protection in and out of your vehicle. ScanX can be operated while sitting in the open case.

Dimensions: 22"L x 22"W x 20" H, PN: D5995

HARD CARRY CASE

Rugged portability is built into this case. Just remove the cover and the ScanX operates while sitting in the lower portion of the case. The outer shell is extremely durable; ABS plastic laminate with ¼" thick plywood base. Thick, sturdy foam secures the ScanX within the case. Three-inch tilt-wheels and telescoping handle help you get this system to wherever you need to work. Recessed twist latches stay securely locked in place. Add your own key-lock or padlock for additional security. Case meets requirements for air travel.

Dimensions: 23"L x 19.8"W x 19.8"H, PN: D5990

WALL MOUNTS

Organization and efficiency is key with our innovative mounting options. Wall Mounts offer a fantastic space-saving solution for your practice. Securely support your valuable investment while freeing up counter space. Choose from two models: Computer and ScanX Wall Mount or the Laptop and ScanX Wall Mount.

Wall Mount for Computer and ScanX, PN: D5950-3

Wall Mount for Laptop and ScanX, PN: D5950-2



Hard Carry Case



Soft Carry Case

scanX

ACCESSORIES FOR ALLPRO SCANX SYSTEMS

CLEANING SHEETS

ScanX Cleaning Sheets are a simple and effective preventive maintenance product to help keep the scanner free of dust. Easy-to-use, just run through the scanner as you would an imaging plate.

PN: B2010, Qty. 12 / PN: B2020, Qty. 25

PSP CLEANING WIPES

PSP Cleaning Wipes for all size Phosphor Storage Plates (PSP) and Plate Protectors. The extra-soft, fabric wipes will not scratch or damage while safely removing dust, hair, dirt and smudges from the imaging surface. These single-use wipes are sure to make clean up fast and easy.

PN: B8910, Qty. 50




PSP Cleaning Wipes

IMAGING PLATES, PROTECTORS AND CASSETTES

ALLPRO Imaging sells a variety of durable, reusable Phosphor Imaging Plates, and the accompanying Plate Protectors and Cassettes in all body and dental sizes: 8" x 10", 10" x 12", 14" x 17" and all intraoral dental sizes, 0 through 4.



Imaging Plates
(pictured in Plate Protectors)



Imaging Cassettes



Scan QR code to learn more about ScanX Accessories

www.allproimaging.com

ALLPRO Imaging (a division of Air Techniques, Inc.)

1295 Walt Whitman Road | Melville, NY 11747-3062

Phone 888.862.4050 | Fax 888.247.8481

9922-310 Rev D

**ALLPRO
IMAGING**
equipped for life™